

Sure Quality is on  
**Sure-endo**



# **SURE-Paste**<sup>®</sup> **SURE-Prep**<sup>®</sup>

Water Soluble Calcium-Hydroxide Paste and  
BLUE EDTA(15%) Preparation Cream from  
**Sure-endo**



**Gutta Percha Points**  
**Absorbent Paper Points**

Hand-Rolled Gutta Percha points and Absorbent Paper Points for Root Canal from Sure-endo

**Just try it**  
feel the differences

Class | Size | Honey  
Control | Size | Variety

**SURE-Cord**<sup>plus</sup>  
Knitted Retraction Cord

Retraction&Hemostatic cord from Sure-endo

**Just try it**  
feel the differences

Class | Spring | Type

**Sure-Seal Root**<sup>™</sup>  
Syringe type & Putty type

Bioceramic Sealing and Filling Material from Sure-endo

**Just try it**  
feel the differences

Class | Size

**Just try it**  
feel the differences

Stable | High | Easy  
Lubricant | Bubble | Blue



**Sure Dent Corporation**


#809, 52, Sagimakgol-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do, Rep. of Korea  
Tel. +82-31-737-4720 E-mail. suredent@suredent.com Web. www.suredent.com



SD-03 | Rev. 10 / 2019

# SURE-Paste®

Water Soluble Calcium-Hydroxide Paste with Barium Sulfate



SURE-Paste® is water soluble Calcium-Hydroxide paste used as temporary filling material for endodontic procedure from Sure-endo.

SURE-Paste® is highly-flowable premixed paste for cleaning, and can be used easily due to water solubility.



- Antimicrobial effect resulting from high pH 12.
- 'WATER-SOLUBLE' with high flowability provides easy filling or removal.

## Application

- Apexification
- Treatment for infected root canals
- Temporary root canal filling
- Internal & External root resorption
- Pulp capping of exposed pulp and vital pulpotomy

## Advantages

- Very Stable without any solidification
- Outstanding antibacterial characteristics
- Easy cleaning and removal due to water solubility
- Easy to apply (Premixed paste & syringe type)
- Impressive biocompatibility
- Highly flowable
- Disposable mm marked tips to avoid a cross-contamination
- Excellent radiopacity

## Composition

- Calcium-Hydroxide
- Barium Sulfate


## Contents

- Regular Kit  
: 2g Syringe × 1ea  
: Length marked Capillary Tip × 20ea  
- Hygienic Packing
- Special Kit  
: 2g Syringe × 2ea  
: Length marked Capillary Tip × 20ea



# SURE-Prep®

Blue EDTA(15%) Preparation Cream



SURE-Prep® is BLUE EDTA(15%) root canal preparation cream from Sure-endo

SURE-Prep® including Ethylenediamine Tetraacetic Acid(EDTA) is a chelating agent which allows efficient cleaning and easier preparation of the root canal.



- 'BLUE' color makes it easy to identify the removal of all of EDTA before obturation.

## Application

- Root canal preparation with files and reamers
- Facile shaping and cleansing of the root canal
- Great lubricating effect
- Provide a clean and non-smear canal wall

## Advantages

- Removal of inorganic substances by chelation(EDTA-Ca)
- Simple application with lubricating effect
- Provide a clean and non-smear canal wall
- Improving the efficiency in injecting into the root canals by a special lubricant.

## Composition

- EDTA 15%
- Urea Peroxide in PPG and PEG 10%

## Contents

- 6g Sryinge × 3ea
- Disposable Tip × 20ea

